



Malvern Spraytec

Laser diffraction









Spraytec training overview

- Introduction to particle sizing
- Introduction to laser diffraction
 - » The diffraction pattern
 - » The Spraytec
 - » Optical models
- Measuring sprays
 - » Multiple scattering
 - » Vignetting
 - » Beam steering
 - » Low signal
- Spraytec measurement process – at the instrument
 - » SOP wizard

Introduction to particle sizing




Particles come in many different shapes (as well as sizes)



- How do we describe the size of these particles?


Concept of equivalent spherical diameters

- Different particle sizing techniques report different equivalent spherical diameters
 - » Dependent on the physical property that is measured




Equivalent sphere theory


- Laser diffraction describes a 3D shape as the diameter of a sphere that would exhibit the measured effect



Number versus Volume





Result = Volume weighted





Introduction to laser diffraction



Dependence of diffraction pattern on particle size




Large particle



Small particles

A basic laser diffraction system...



Spraytec system detail



Spraytec Lens Ranges





- Wide Dynamic Range
 - » 300mm lens: 0.1 – 900 microns
 - » 750mm lens: 2-2000 microns (5 – 1600 microns for Dv50)
- Large working distance
 - » 150mm at 0.5 microns
 - » Over 1000mm at 5 microns
- Small Surface Area
 - » Minimises exposure
- Purge system available for large volume sprays




Spraytec Data Capture




- The collected diffraction pattern is converted to a Volume Particle Size Distribution (equivalent sphere) following the Mie theory (or Fraunhofer)




Spraytec dynamic particle size analysis



How the size distribution is calculated?




Fraunhofer approximation



Scattering data analysis


Mie Theory



“..... the Mie theory offers the best general solution.”


ISO 13320-1

Scattering Data Analysis Using different optical models



Blue Fraunhoffer
Red Mie


Scattering Data Analysis Using different optical models




Measuring sprays



Particle size by Spraytec




Spraytec Measurement Process




Alignment + Background

Spraytec Measurement Process



Measurement


Spraytec Measurement Process



Results!

- Measurement manager ensures each measurement step occurs in the correct order, improving result reliability

Spraytec Result Interpretation




Spraytec Result Interpretation




Size History

- Frame-by-frame event analysis
- Drag-and-drop result selection
- Zoom facility to examine changes over different time periods
- Over plotting of results at different time points
- Display of measurement triggers and data groups



Background


- Background signal should decrease across increasing detector numbers
- As the light energy from the laser decreases with increasing angle



Background quality

■ Droplet deposition on the windows

- » Unfortunately sprays are messy, and not always easy to control so droplets can deposit on the windows.




Preventing deposition on the windows

- An air purge can be used for particularly messy sprays.
- Otherwise, a background can also be measured after each measurement.
- An acceptable % variation of the background level can then be set.
- If this % change is exceeded then a warning to clean the windows will be displayed to the user.




Background warnings

- The background warnings can be set up in the SOP on the Measurement setup page



Ambient light problems


- Background measurements can also be contaminated by ambient light
- This is seen as fluctuations on the outer channels and can cause problems when signal levels are low.
- Check location of the instrument: nearby window? Other light source?



Spraytec Multiple Scattering




- Multiple scattering causes the measured scattering to move to higher angles, yielding too fine particle size.
- In the Spraytec software it is possible to correct for this.



Spraytec Multiple Scattering


- Example : higher concentration measured and multiple scatter occur
- -> extra tail of small particles (red)



Spraytec Multiple Scattering



- In the Spraytec software it is possible to correct for this.




Vignetting

- High angle scattering data can be lost if the spray is positioned too far away from the detector lens.
- The maximum distance at which a spray can be measured is known as the working distance.
- This is dependent on the particle size.
 - » 150mm for 0.5 micron particles
 - » Over 1000mm for 5 micron particles
- Loss of high angle scattering can cause the measured particle size to increase, or the distribution to narrow.

Spraytec Lens Operation




Vignetting



Working Distance depends on the particle size (scattering angle) and the lens diameter


Vignetting: Inspecting the scattering data

- Scattering data from a spray affected by vignetting
- The data can be reanalysed without the affected data channels




Vignetting

- This graph shows the particle size distribution for the same spray measured at increasing distance from the detector.
- Vignetting begins to occur at distances over 15cm.




Spray position and vignetting



Vignetting – remove data channels


Removing data channels

- Then the data can be reanalysed using the reduced detector range
- Edit records and then alter the detector range (can also be set up in a SOP)



Vignetting: Remove data higher detectors


- Remove data higher detectors
 - » Measurement range limited




Beam steering: RI change and beam steering

f.e. Sprays contain volatile components used as propellants.

Unexpected peak large particles.



Beam Steering




Beam Steering




Removing data channels

- Then the data can be reanalysed using the reduced detector range
- Edit records and then alter the detector range (can also be set up in a SOP)




Beam Steering



Velocity Bias


What happens if the extraction system is not effective:

- An extraction system is needed to ensure that all droplets travel through the laser beam at the same velocity.
- Without extraction the fine particles decelerate more quickly and are over sampled.
- An inefficient extraction system can cause a reduction in the measured particle size.




Velocity Bias

- Effect on the measured size
- Efficiency of the extractor is reduced with distance and the particle size decreases.




Velocity Bias

- Effect on the profile
- If the extraction is not sufficient droplets may also fall back through the beam, affecting the stable phase of the transmission profile




Measurement at low and high concentrations

- Seeing red?
- If you have measured a spray and you see no results, only red bars, you may have to adjust the warning levels.




Measurement at low and high concentrations

- Warning and Error levels.
- These can be altered in case you have a spray with a particularly low or high concentration




Measurement at low and high concentrations

- Applying the new levels
- By setting the warning and error levels to zero, you will be able to see the data
- From which you can set a more appropriate error level to be used in future measurements.



Measurement at low and high concentrations

- Low signal levels
- If you are measuring at low signal levels, then it is also important to check the signal threshold level.
- Any data below this threshold will not be analysed this will affect the quality of the fit.



METHOD DEFINITION AT THE INSTRUMENT






- SOP Wizard
 - » Similar to the Mastersizer 2000
- Guides users through the process of method development
- Automates the measurement process
- Template SOPs provided for different spray types.

Spraytec SOP Settings

SOP Type



■ Hardware configuration

- » Includes accessories

► Basic measurement settings:


- Continuous mode
- Rapid mode
- Timed mode

► Lens Range

- 300mm lens
- 750mm lens

Spraytec SOP Settings

SOP Type



■ Hardware configuration

- » Includes accessories

► Basic measurement settings:

- Continuous mode
- Rapid mode
- Timed mode

► Lens Range

- 300mm lens
- 750mm lens

Spraytec SOP Settings Labels and Instructions



- Similar functionality to the Mastersizer 2000
- Provides default sample details
- Prompts the user to carry out specific actions as part of running an SOP

Labels

Sample name:
My Spray

General comments:
Direct spray towards somebody else

Source type: Works Source name: Nebs "r" Us Lock source

Lot type: Run Lot value: Away Quickly Lock lot Auto-increment

Show this page when the SOP runs, so that the user may modify and add notes.

< Back Next > Cancel Help


Instructions

Display the following instructions before the measurement begins.
Wear raincoat to protect Armani Suit

Display the following instructions after the measurement ends.
Smile - you know it was all worth it.

< Back Next > Cancel Help

Actuator Settings



Actuator type


- » InnovaSystems or Proveris

Spraytec NSS settings

- Distance and angle
- Requirements for reference setting

Spray profile

- Aids the concentration calculation
- Leave as default unless Cv is important to you...



Cell Angle Setting

- » Set vertical, 45 degrees or horizontal
- » This is not verified during the measurement


Extraction Settings

- Controls the primary pump
- Switch on: At alignment
- Switch off: At end of SOP

Spraytec SOP Settings



Measurement Setup




- Background duration
- Background drift warning
 - Enables background check after measurement
 - Determine %age draft across the selected range of detectors
- Alarm settings
- Enable Output trigger
 - Synchronise other devices with the measurement start
- Enable Inspection Stage
 - Pauses the measurement after the background measurement

Spraytec SOP Settings



Measurement Setup: Alarms



Set Default values

- Will ensure measures occur apart from if a significant error occurs

Tests for the optical background

- Minimum transmission signal
- Average scattering
- Peak scattering signal


Tests for the electronic background

- Average scattering
- Peak scattering signal

Background expiry time

Spraytec Measurement Triggering

Rapid mode



- **Triggering based on:**

- » Manual mouse click
- » Transmission level
- » Scattering level
- » External Signal
 - TTL trigger
 - Switch trigger


- **Can collect data from before the trigger event**

- » Max buffer size of 3 sec

- **External trigger can also be used to group data**

- » Useful when multiple events occur within one measurement

Rapid measurement Advanced Settings



Filters

- **Transmission:** Removes all data above a specified transmission level
- **Signal:** Removes all data below a specified average signal level


Can put data into groups based on an external trigger

- Only used for rapid, repeating, pulsed spray events

Spraytec SOP Settings



Analysis Options



Particle Properties

- Refractive Index
- Absorption

Dispersant Properties


- Refractive Index
- Normally set to Air

Data handling options

Multiple scattering analysis

- Switched on as default
- Should be disabled for nasal sprays

Analysis Options: Data Handling



Detector range for the analysis

- Allows data to be removed
- Can apply first ring kill using an automatic algorithm

As default, the extinction analysis is turned off

- Requires an accurate path length to work...


Scattering Threshold

- Data required before a detector channel is included in the analysis
- Normally set at 0 or 1

Spraytec SOPs Analysis Options




Data Handling : Spray Profile



- Input path length or spray profile information
 - » Important for concentration calculations and extinction analysis

Spraytec SOPs Analysis Options

Data Handling : Analysis Alarms



Set Default values

- Will ensure measures occur unless a significant error occurs

Minimum transmission alarm

- Stops size calculation if the spray is too concentrated


Minimum scattering alarm

- Stops size calculations if the spray is too dilute

Spraytec SOPs Analysis Options


Data Handling : Analysis Alarms

Sample : Placebo nasal dpi
 Event: 1 +0.0010 (s)



Title	Legend	I	Correction Equation
Trans (%)	—	—	—
Dv(10) (μm)	—	—	—
Dv(50) (μm)	—	—	—
Dv(90) (μm)	—	—	—

Additional Parameters



Use 2-parameter curve fit

- Rosin Rammler model
- Log-Normal model
- Need to specify percentile range for the curve fit


Change User sizes

- Size range
- Linear or log spacing

Advanced Options

- Allows the default size range of the analysis to be changed

Spraytec SOPs Derived Parameters



► List of parameters which will be calculated as part of a measurement

- Can define a maximum of 16


► Add new parameters here....

► Can set an averaging period for each parameter

- Important if you decide to display the average or StDev on the size history display.

Spraytec SOPs Rapid Measurement

Averaging



■ Data-based average

- » Calculated by averaging the **scattering data** over a specific time period. A new size distribution is then calculated for this **average data set**

■ Concentration weighting


- » When enabled, corrects for changes in the particle concentration when calculating an average

■ Result-based average

- » Calculated by averaging the **size distribution size band data** over a specific time period, producing an average size distribution

Spraytec SOPs Rapid Measurement

Time-based Averaging




- Averages the data over fixed time periods
- Time period specified by
 - » External trigger events
 - » Specified time interval

- Produces a new Particle Size History containing the average results from each time period

Spraytec SOPs Rapid Measurement


Phase Averaging



- Averages the data across a maximum of 3 separate “phases” with an event
 - » Start / Formation Phase
 - » Middle / Fully Developed Phase
 - » End / Dissipation Phase
- First phase starts at the beginning of the event
- Produces average size distributions for each phase and a new PSH if required
- Meets the FDA's guidance for nasal spray testing


Nasal sprays

Atomisation dynamics



Nasal sprays


Atomisation dynamics



Nasal sprays




Delivered particle size



Phase	D _{v10} / μm	D _{v50} / μm	D _{v90} / μm
Formation Phase	34.54	67.93	690.13
Fully Developed Phase	19.33	37.2	58.31
Dissipation Phase	31.83	45.27	75.36


Spraytec SOPs Rapid Measurement

Ensemble Averaging



- Use to calculate an average event
- Averages the data for the same time point in each event
 - » Start and stop time defined by the trigger position
- Produces an average Particle Size History file

Spraytec SOPs Wizard Complete!





Spraytec Inhalation cell



PHARMACEUTICAL AEROSOLS

When can we use Laser Diffraction?




- Nebulizers and soft mist inhalers
- Dry powder inhalers (DPIs)
- Meter-dosed inhalers (MDIs)
- Nasal sprays
- What might we need?
 - *Inhalation cell*
 - Pump
 - Flow meter
 - Critical flow controller
 - *Nasal spray support*
 - Actuator

Pharmaceutical Aerosols

Why use Laser Diffraction?

- Access to Dynamic Information
- Formulation Screening Capabilities
- Measurements under “real” conditions



STP2300 Spraytec Inhalation cell



- For the measurement of
 - » Metered-dose Inhalers
 - » Dry-powder inhalers
 - » Nebulisers
- Enables measurements in an enclosed cell at a controlled flow rate
- Can be mounted
 - » Vertically – measurement in the throat
 - » Horizontally – measurement in the mouth
 - » 45 degrees – nebulisers
- Compatible with standard inhaler testing equipment (e.g. USP throat, cascade impactor)



summary regulatory requirements for DPI particle size

- Each device has a different resistance to air flow.
- Flow rate is defined to give a certain pressure drop of across the device
 - » USP requires 4kPa (with max flow rate of 100l/m)
- The inhalation volume is also defined:
 - » USP : 4L
 - » FDA/Eur. Ph. : 2L
- Critical flow control is required
- Pre separator required for cascade impaction

STP2310: USP Induction Port




- USP Throat from Copley Scientific
- Generally used with
 - » Metered-dosed inhalers
 - » Dry-powder inhalers



Andersen Cascade Impactor ACI



Andersen Cascade Impactor ACI



STP2311-(220/110) HCP5 High Capacity Pump



- From Copley Scientific
- Designed for testing inhalers
- Regulated vacuum from 0l/min to 250l/min
 - » Dependent on inhaler device resistance
- Unregulated output of 416l/min
 - » For use with flow controllers(TPK2000) for DPI testing
- Requires a switching box to be synchronised with the Spraytec
- This pump is designed to be used with devices which deliver particles over short time periods – not minutes



STP2329 TPK2100 critical flow controller



- Flow control valve
 - » Automatically opened using the Spraytec TTL trigger
 - » Control the flow use to activate DPI
 - Desired pressure drop is achieved instantaneously at the start of the measurement
 - Monitors pressure drop across the flow valve to ensure sonic flow is achieved
- Measures pressure drop across an inhaler
 - » Requires STP2305 – dose collection device
 - » 0kPa to 15kPa
- Can log flow data from DFM2000 (STP2314)
- STP2328 Copley Scientific Calibration Kit for the TPK2100 critical flow controller



STP2314 DFM2000 Flow meter




- Measures flow rate produced by the HCP5 pump (STP2311)
- Measurement range
 - » 0l/min to 290l/min (0.1l/min resolution)
- Interfaces with TPK2000 critical flow controller
- Includes an adapter to attach the flow meter to the UPS induction port



STP2305 Dose Collection Device



- From Copley Scientific
- Collects the dose emitted from an inhaler prior to a vacuum pump
 - » Prevents damage to the pump
- Required during DPI testing with the TPK2000 critical flow controller
 - » Measures pressure drop across the device
- Includes 100 glass filters
 - » And replacement seals




STP2306 Mouth Piece adapter

- from Copley Scientific
- Connects inhaler to
 - » USP Throat (STP2310)
 - » Dose collection device (STP2305)
 - » Or directly to the inhalation cell
- Requires purchase of moulding tool STP2307
 - » Example inhaler must be supplied with the order



Metered-dose inhaler set-up




Dry-powder inhaler set-up



Input:

USP induction port (STP2310) &
Mouthpiece adaptor (STP2306 &2307)

Flow meter
DFM2000
STP2314



Spraytec and inhalation cell:
STP5315

Trigger output
Cable:

Particle collection:
Cascade impactor
or Dose collection
unit (STP2305)

Digital output
Cable:

Critical flow controller
TPK2000 (STP2313)

Switch box


Pump:
HCP5(STP2311)

Nebuliser set-up



Input:
Mouthpiece adaptor (STP2306 &2307)
T-piece (customer)

Flow meter
DFM2000
STP2314



Spraytec and inhalation cell:
STP5315 in horizontal position

Digital output
Cable

Switch box
STP2400

Particle collection:
Filter (STP2309) or
collection unit
(STP2305)




Pump:
HCP5(STP2311)

Dry Powder Inhalers

Principle of operation


- Delivers micronised drug (+ excipient)
 - » Drug: 1 -10 microns
 - » Excipient: 60 - 200 microns
- ▶ Drug Delivery Process:
 1. Powder is metered prior to delivery
 2. Powder is entrained when patient inhales
 3. Shear is applied to disperse agglomerates
 4. Aerosol is delivered to the patient



Pharmaceutical Aerosols Why?

Inhaled drug formulations


- Throat Deposition
 - » Impaction ($>5\mu\text{m}$)
- Upper Airways
 - » Sedimentation (1-5 μm)
- Lung Deposition
 - » Diffusion (<2 μm)



Clive Washington: Particle size analysis in pharmaceutics and other industries:
 Theory and Practice p179
 Ellis Horwood, 1992 ISBN: 0 13 651613 0


Particle size flow rate dependence

Unprocessed lactose powder




Particle size flow rate dependence

Grindomix powder




Particle size flow rate dependence

Mechofused powder



Changes to process


Changes aerosolisation behaviour...



Changes to process

Changes in aerosolisation behaviour...


Agglomerates



High Shear Blender




Agglomerates are
coated



+
Mg Stearate


Mechanofusion



Each grain is
dispersed and coated

Cv time dependence


Mechofused powder



Dry Powder Inhaler




Comparing Formulations



Dry Powder Inhaler




Comparing Formulations




Metered Dosed Inhalers

- Canister containing drug substance in suspension or solution
- Solvents or surfactants are added to aid solubility or dispersion.
- Solution is metered by the valve.
- Aerosol is produced as the solution is forced through the metering valve by the propellant.



Example of measurement

- Likely to encounter beam steering due to the propellant.





summary of regulatory information

Metered Dose Inhalers


- Flow rate of 28.3l/m for use with ACI
- Pre-separator is not required
- Critical flow apparatus is not required

Nebuliser Characterisation



General Overview


- Drug dissolved in aqueous solution
 - Liquid shearing produces droplets
 - » Primary droplets are re-circulated
 - » Satellite droplets are inhaled
 - Important parameters include:
 - » Jet and baffle design
 - » Surface tension and viscosity of drug solution
 - » Gas pressure, density and flow rate
- Good correlation between laser diffraction and impaction results
- Avoids evaporation problems seen with inertial impactors
 - Rapid characterization



Particle size distribution


Beam steering can be observed from nebulisers

- » Due to the evaporation of very fine droplets
- » Smaller effect than that observed in MDI's




Nebuliser Characterisation

12 l/min driving gas flow rate



Nebuliser Characterisation Measurement During Inhalation





Spraytec Nasal Spray



Spraytec Nasal Spray Measurements

- Spraytec NSS Nasal Spray Support
- Vertical Position Control
 - » 0 – 100mm
- Control of angle
 - » 0 – 45 degrees
 - » 5 degree steps
- Position feedback to software

Regulatory perspective for nasal spray testing

- Device and formulation interaction are critical
 - » Nasal sprays are combination products
 - » Need to ensure reproducible delivery
- Droplet size analysis is used to assess
 - » Quality: Product performance
 - » Safety: Lung penetration
 - » Efficacy: Bioavailability and Bioequivalence Studies for Deposition Alveoli and Nasal Sprays for Local Action



U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)

April 2003

Nasal sprays: regulatory perspective



- Device and formulation interaction are critical
 - » Nasal sprays are combination products
 - » Need to ensure reproducible delivery
- Droplet size analysis is used to assess
 - » Quality: Product performance
 - » Safety: Lung penetration
 - » Efficacy: Deposition location




STP2213 Spraytec NSS nasal spray support

- References a Proveris Vereo or Proveris Solo nasal spray actuator relative to the measurement zone
 - » Feedback of distance between nozzle and laser beam (0-10cm range)
 - » And angle of actuation (vertical – 45 degree range)
- Measurement distance and angle can be specified in the Spraytec SOP
- Only the Solo actuator (velocity controlled) can be supplied by Malvern



STP2214 Spraytec NSS nasal spray support

- References an Innova systems NSP3000 Mighty Runt nasal spray actuator relative to the measurement zone
 - » Feedback of distance between nozzle and laser beam (0-10cm range)
 - » And angle of actuation (vertical – 45 degree range)
- Measurement distance and angle can be specified in the Spraytec SOP
- The Mighty Runt (force controlled) must be purchased directly from Innova Systems



STP2325 Solo for Spraytec nasal spray actuator

- Controls velocity of actuation
 - » Characterises the bottle size
 - Seats bottle before actuation
 - » User inputs
 - Maximum velocity

For actuation and return

- Hold time
 - Delay between actuations
 - Trigger mechanism
- » Outputs
 - Force vs time profile
 - Position vs time profile
 - Option to record shot weight



Measurements specified by the FDA


Nasal Sprays

- Assessment of the **spray dynamics**
 - » Provide “size history” of the spray event
 - » Allows Formation, Stable (Fully Developed) and Dissipation phases to be observed
- Provide following parameters for the **Fully Developed Phase**:
 - » Dv10
 - » Dv50 (Median)
 - » Dv90
 - » Span ((Dv90-Dv10)/Dv50)
 - » (%<10 microns)
- Measurements required for the **B and E life stages**
- **Distance dependence** need to be studied
 - » Distance of between 2cm and 7cm
 - » Separation of at least 3cm required.
- **3 repeats** required for each measurement to assess reproducibility
- Use of **automated actuation stations** is recommended

Nasal sprays




Formulation viscosity



Nasal sprays

Formulation viscosity





Nasal spray set-up

- Spraytec, short bench, 300mm lens (STP5311)
- Nasal spray support, NSs (STP2213 or 2214)
- Actuator
 - » Proveris Solo (supplied by Malvern, STP2325)
 - » Proveris Vereo
 - » Innova systems
- Extractor (if none available in the lab)
 - » STP2499 (110v) & STP2500(220v)



Spraytec extraction system




STP2499 (110v) & STP2500(220v) Extraction system

- For the extraction of low volume sprays

- » Domestic aerosols,
- » nasal sprays
- » and personal products.


- Includes

- » spray collector assembly,
- » tubing,
- » HEPA filter unit
- » switching box to synchronise with the Spraytec





Spraytec fuel injection application



Fuel injection





Pulsed sprays: Fuel injection systems


- Particle size is related to:
 - » Combustion efficiency
 - » Production of harmful emissions
- Small droplet size:
 - » Increases available surface area
 - Faster evaporation
 - » Decreases spray penetration
 - Affects mixing process
- Important parameters include:
 - » Injection pressure
 - » Injector type / geometry
 - » Valve timing / operating speed

Fuel injection experiment

- Rapid
 - » 10kHz (needs software key)
 - » Minimum scattering alarm: 5
- Short events
 - » 10ms
- Multiple events
 - » 5-20
- Water or fuel
- Possible beam steering from evaporation

Fuel injection


Single injection event



Fuel injection


Affect of fuel pressure

Cumulative volume %




Fuel injection

Assessing pulse variation



Fuel injection

Effect of pressure on variability






Spraytec agrochemical sprays application

Continuous sprays

Agricultural sprays

- Sprays enable efficient coverage of large areas.
- Why is droplet size important
 - » Targeted absorption
 - Finer droplets deposit on leave
 - Larger droplets penetrate the canopy
 - » Reduce risk of run off
 - Associated with large particles
 - » Reduce the risk of spray drift
 - Associated with small particle
- Nozzles are characterised to match target and environment






Agrochemical spray experiment

- 750mm lens + long bench
 - » Can use 300mm lens for demo if necessary
- Continuous
 - » Max 1 hour at 1 Hz
- Water or formulation
- Possible beam steering from air currents with long benches
- Protect alignment by using shrouds to prevent liquid settling on the X95


Agrochemical sprays

BCPC boundary nozzles




Agrochemical sprays

Nozzle traverse




Agrochemical sprays

Affect of pressure



Classification of agrochemical nozzles

- Size boundaries defined by BCPC nozzles
- Each nozzle has been measured over a range of pressures
- Nozzles can then be classified according to their size over standard pressure ranges





Spraytec household and domestic sprays application



Personal and domestic products

- Very large market for spray systems
- Products include:
 - » Perfumes
 - » Deodorants
 - » Hair Sprays
 - » Furniture Polish
 - » Air Fresheners
 - » Cleaning products
 - » Pest control
- Particle size important in order to:
 - » Ensure even surface coverage
 - » Provide rapid evaporation / drying
 - » Reduce “wetness”
 - » Prevent lung deposition (no particles < 10µm in size)

Aerosol experiments


■ Continuous

- » Acetone or standard opaque particle
- » Beam steering from propellant

■ Rapid


- 100Hz – 2.5kHz
- » Medium-long events
 - 500ms – 5000ms
- » Multiple events
 - 3 – 5
- » Water
- » Possible beam steering from evaporation

Hairspray formulations



Personal products

Hairspray formulations



Spraytec wet dispersion



Wet dispersions

